



TIMS DICOM[®] System Integrating the Healthcare Enterprise (IHE) Integration Statement

**Revision 1.0
October 13, 2006**

Copyright © 2006 by Foresight Imaging[®]. All rights reserved. All specifications are subject to change without notice.

This IHE Integration Statement describes the intended conformance of the TIMS DICOM System with the IHE Technical Framework. It identifies the specific IHE capabilities this product is designed to support in terms of the key concepts of IHE: Actors and Integration Profiles.

| IHE Integration Statement | | | |
|---|--|---|-------------|
| Vendor | Product Name | Version | Date |
| Foresight Imaging | TIMS DICOM System | Version 2.0.4042 | 10/13/2006 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Scheduled Workflow | Acquisition Modality (Query Modality Worklist completed, all other transactions under development.) | None | |
| | Image Display | None | |
| Internet Address for Vendor's IHE Information: www.tims.com/ihe.htm | | | |
| Links to Standards Conformance Statements for the Implementation | | | |
| DICOM | www.tims.com/dicom-conformance-statement.htm | | |
| Links to General Information on IHE | | | |
| North America: www.rsna.org/IHE | Europe: www.ihe-europe.org | Japan: www.jira-net.or.jp/ihe-j | |

Foresight Imaging
 70 Industrial Ave.
 Lowell, MA 01852
 Tel: 978-458-4624
 Fax: 978-458-5488
 Web: www.tims.com
 Email: info@tims.com